

## **Certificate in relation to air transportation by passengers with walking disabilities**

We hereby certify that our product

### **e-motion M25 (incl. DuoDrive)**

including the lithium-ion batteries fulfil the requirements of

### **IATA chapter 2.3.2.4**

and are therefore suitable for transportation by passengers according  
The pilot in command must be informed about the location of the mobility aid.

The built-in lithium-ion batteries (**150 Wh**) of the e-motion M25, which are not intended for removal by the user, and the DuoDrive control unit (**< 10 Wh**) fulfil the requirements of the tests listed in the UN Manual of Tests and Criteria" Part III subsection 38.3.

The e-motion M25 can be set into a dedicated "flight mode". In this mode the system is completely deactivated. While in this mode the e-motion M25 can only be switched on by a coded procedure. Unintended switching on by luggage is not possible. Furthermore, the Bluetooth device of the system is completely switched off. The DuoDrive control unit must be taken into the cabin with the hand luggage.

**We strongly recommend contacting the airline sufficiently in advance to request approval for transport.**

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31<sup>st</sup> December 2024**.

Albstadt, 26<sup>th</sup> of October 2023

Alber GmbH



Joost Beltman  
Managing Director