

Certificate

We confirm, that the nickel metal hydride batteries used in our product

e-motion (M10, M12)

are classified as

UN Nr. 2800 „Batteries, gel, nonspillable, electric storage“

according to IATA (International Air Transportation Association).

According to special provisions A-67 of IATA Dangerous Goods Regulations, they are considered as non-dangerous goods on transportation. No electrolyte will flow from cracked case at a temperature of 55°C.

The terminals have to be protected against short circuit by the delivered covers!

They are also considered as „nonspillable battery“ by DOT (U.S. Department of Transportation).

Albstadt, 2005-01-03

Ulrich Alber GmbH + Co. KG



Ralf Ledda
Managing Director