

ZERTIFIKAT ♦ CERTIFICATE ♦ CERTIFICADO ♦ CERTIFICAT ♦ CERTIFIKAT ♦ 認証証書 ♦ CERTIFICATE ♦ CERTIFIKAT ♦ CERTIFICADO ♦ CERTIFICAT



Product Service

CERTIFICATE

No. Z1 05 10 12941 101

Holder of Certificate: Invacare Deutschland GmbH

Kleiststrasse 49
32457 Porta Westfalica
GERMANY

Production Facility(ies): 12941

Certification Mark:



Product: Wheelchair for handicapped

Model(s): Electrically driven wheelchair type: G50

Parameters:	class of use:	C
	total height:	1065 mm
	total length:	1240 mm
	total width:	700 mm
	max. load:	150 kg
	max. safe slope:	10°
	max. obstacle height:	100 mm
	driving motor:	Invacare GmbH, type: DG2022
	servomotor	
	(steering support):	Invacare GmbH, type: S141054
	electronic:	Dynamic Controls Ltd., type: DX
	batteries:	MK Battery, 2x (12 V, 83 Ah (5hr), type: M24 SLD GFT); 2x (12 V, 80,5 Ah (5hr), type: M31SLDG)
	battery charger:	CTE Corporation, model: BAT-GCO812, 24 V DC, 8 A; POWERCHARGE / GI-TEAMTEC, type: GBC04, 24 V DC, 12 A

Tested according to: EN 12184:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. See also notes overleaf.

Test report no.: 70109263

Date, 2005-11-09

Page 1 of 1



186808



Product Service

CERTIFICATE

No. Z1 05 10 12941 100

Holder of Certificate: Invacare Deutschland GmbHKleiststrasse 49
32457 Porta Westfalica
GERMANY**Production Facility(ies):** 12941**Certification Mark:****Product:** Wheelchair for handicapped**Model(s):** Electrically driven wheelchair type: G50

Parameters:	class of use:	B
	total height:	1085 mm
	total length:	1240 mm
	total width:	700 mm
	max. load:	150 kg
	max. safe slope:	10°
	max. obstacle height:	100 mm
	driving motor:	Invacare GmbH, type: DG2022
	servomotor	
	(steering support):	Invacare GmbH, type: S141054
	electronic:	Dynamic Controls Ltd., type: DX
	batteries:	MK Battery, 2x (12 V, 63 Ah (5hr), type: M24 SLD GFT); 2x (12 V, 80,5 Ah (5hr), type: M31SLDG)
	battery charger:	CTE Corporation, model: BAT-GC0812, 24 V DC, 8 A; POWERCHARGE / GI-TEAMTEC, type: GBC04, 24 V DC, 12 A

Tested according to: EN 12184:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. See also notes overleaf.

Test report no.: 70109263**Date,** 2005-11-07

Page 1 of 1



186806