



# INVACARE® TOP END® PRO™ TENNIS WHEELCHAIR



For new players, adjustability is key. The **Invacare® Top End® Pro™ Tennis Wheelchair** has just that, featuring adjustable front and rear seat heights, center-of-gravity, footrest angle/height, backrest angle/height and anti-tip. Plus, the Invacare Top End Pro Tennis Wheelchair is available in two frame height choices: short or tall. The short frame is ideal for children and the tall frame is for adults, which makes this chair a great choice for programs that have multiple users. Adjusts in minutes, super easy to use, lightweight and durable. This chair is a winner.

## Features

Adjustable front and rear seat heights  
Anti-tip/5th wheel: Adjustable height  
Fixed back  
Adjustable back height  
Adjustable tension back upholstery  
Adjustable Center-of-Gravity  
Fixed camber 15°, 18° or 20°  
High performance Sun® wheels:  
24", 25", 26"

Anodized pushrims  
High pressure clincher tires  
Smart fork caster system  
72mm casters  
Adjustable platform footrest  
Powder coat paint  
250 lb. weight capacity

## Model #: SPTN



1. Front seat height and center-of-gravity adjustments
2. Rear seat height adjustment
3. Quick release lever
4. Single swivel fixed fifth wheel
5. Precision fork caster system
6. Adjustable footrest system



**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

### Top End

www.topendwheelchair.com  
Sports and Recreation Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2014

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and ®.  
All trademarks are owned  
by or licensed to Invacare  
Corporation unless  
otherwise noted. Spinerger  
is a registered trademark of  
Spinerger Inc. Sun and Equalizer  
are registered trademarks of  
Sun Components. Kenda is  
a registered trademark of  
Kenda Rubber Industrial  
Co., Ltd. Specifications are  
subject to change without  
notification.  
Form No. 12-507  
Rev. 12/14  
141157

## Options

Clamp-on fixed aluminum side guards  
Spinerger® SLX wheels  
Sun® Equalizer™ wheels  
12 tab pushrims  
Kenda® tires  
Click straps  
O-ring positioning belt

Spoke guards  
Frame protectors  
Alignment gauge  
Telescoping anti-tip bar  
Cushion  
Folding Back  
Fold-in side guards

