



*Invacare*<sup>®</sup> | **Comet<sup>®</sup> Alpine+**

Challenge the outdoors with power,  
control and stability



For those who enjoy adventures across terrains too challenging for a traditional scooter, the **Comet<sup>+</sup> Alpine<sup>+</sup>** offers the perfect solution. Its robust and powerful driving characteristics gives you the freedom to discover the great outdoors with minimal effort in continuous comfort and safety.

As part of the **Comet** family, the **Comet<sup>+</sup> Alpine<sup>+</sup>** comes packed with all of its well-known safety features and customisation capabilities, offering a reliable, powerful and stylish scooter.

- Powerful and robust
- Excellent outdoor capabilities
- Safety at its core

## Powerful and robust

Whether you're tackling uneven terrain or steep and challenging hills, the **Comet<sup>+</sup> Alpine<sup>+</sup>** will offer excellent maneuverability and power without compromising safety and style. Thanks to its powerful motors and superior suspension system, the **Comet<sup>+</sup> Alpine<sup>+</sup>** lets you challenge the great outdoors with minimal effort and complete control.

## Features and options



### Seat suspension module

Offers additional comfort by adding an optional seat suspension.



### Seat riser

Electronically increase the seat height to achieve a preferred driving height.



### Off-board charging

An optional off board-charger allows you to conveniently charge the batteries away from the scooter.

## Excellent outdoor capabilities

When faced with wet or uneven ground, you can rest assured that vital components will be protected from water and dirt thanks to the in-built anti splash guards. Feel stable and at ease when driving at high speeds with directional stability to keep you moving in a straight line and automatic speed reduction to slow you down when driving round corners.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Backrest width	Seat depth	Seat to floor height	Backrest height	Driving unit width	Total length	Empty weight	The heaviest part	
<b>Comet<sup>+</sup> Alpine<sup>+</sup></b>	510 mm	Comfort and Deluxe: 470 mm Premium: 460 mm	Min. 440 mm Max. 510 mm	Comfort: 500 mm Deluxe: 490 mm Premium: 630 mm	665 - 685 mm	1450 mm	143 kg	approx. 63 kg	
	Max. user weight	Battery capacity	Turning diameter	Max. climbable obstacle height	Ground clearance	Max. safe slope	Max. driving distance	Max. driving speed	
<b>Comet<sup>+</sup> Alpine<sup>+</sup></b>	160 kg	2 x 12 V/75 Ah (C20) AGM	2750 mm	100 mm	100 mm	10°	49 km	10 km/h	
	Motor capacity	Colour of shrouds			Upholstery				
<b>Comet<sup>+</sup> Alpine<sup>+</sup></b>	10km/h: 650 W / 1800 W <sup>1</sup>	Marble White	Sapphire Blue	Jade Green	Luxor Beige	Crystal Copper	Black Vinyl		
		Carmine Red	Iconic Chrome	Diamond Silver	Electric Blue	Jasper Red			

1. Peak performance

Tel: +41 (0)22 354 60 10  
 Fax: +41 (0)22 354 60 11  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International Sàrl  
 Route de Cité-Ouest 2 - 1196 Gland  
 Switzerland

Trusted Solutions,  
 Passionate People

