



Yes, you can.®

Januari 2021

Invacare nv
Autobaan 22 - 8210 Loppem - België - www.invacare.be

Invacare® **Top End®**
Force RX Handcycle



Datum :

Referentie:

Dealer:

Gelieve uw keuze aan te geven door een kruisje in het vierkantje

Standaard optie

Optie zonder meerprijs

CODE	Omschrijving	Publieke prijs in €, 6% BTW inbegrepen
BASE MODEL:		
TE10009	Force RX Handcycle	8.000,00
Assembly:		
TE10046	Unassembled	
TE10047	Assembled	
WHEEL SIZES:		
TE50170	Front and Rear; 26"	
TE50172	Front and Rear; 700C	
TE50171	Front, 26" and Rear; 20"	450,00
WHEEL TYPE:		
26" FRONT AND 26" REAR:		
TE50183	Omit Wheels & Axles, Build for 26" Front and 26" Rear	
TE50173	High Performance, 26" Front High Performance, 26" Rear (Clincher)	
TE50174	Corima Aero Carbon, 26" Front Corima Aero Carbon, 26" Rear (Tubular)	3.900,00
700C FRONT AND 700C REAR:		
TE50185	Omit Wheels & Axles, Build for 700C Front and 700C Rear	
TE50181	Corima Aero Carbon, 700C Front Corima Aero Carbon, 700C Rear (Tubular)	3.900,00
26" FRONT AND 20" REAR:		
TE50184	Omit Wheels & Axles, Build for 26" Front and 20" Rear	
TE50176	High Performance, 26" Front High Performance, 20" Rear (Clincher)	310,00
TE50177	C7 Carbon, 26" Front C7 Carbon, 20" Rear (Clincher)	3.250,00
TE50178	Corima Aero Carbon, 26" Front Corima Aero Carbon, 20" Rear (Tubular)	3.900,00
TE50179	Corima Aero Carbon, 26" Front Corima Aero Carbon, 20" Rear (Tubular)	5.200,00
WHEELBASE: (see definitions)		
TE75014	Standard Length	
TE75019	Extended Length	
SEAT WIDTH: (see diagram I)		
TE15011	13" Seat Width	
TE15013	14" Seat Width	
TE15015	15" Seat Width	
TE15017	16" Seat Width	
TE15019	17" Seat Width	
TE15021	18" Seat Width	
CAMBER ANGLE:		
TE45001	0 Degree Camber	
TE45002	3 Degree Camber	



CODE	Omschrijving	Publieke prijs in €, 6% BTW inbegrepen
FORK OPTIONS: (See Definitions)		
FORK STYLE:		
TE65181	Adjustable Height Fixed Length	
TE65180	Adjustable Height Adjustable Length	150,00
LOWER FORK STYLE:		
TE75025	Aerodynamic	
TE75026	Wide Required with Disk Brake	
SHIFTER OPTIONS:		
TE65224	Omit All Shifting Components (see definitions)	
10 SPEED:		
TE65228	SRAM GX Trigger Shifting Package Shifter & brake mounted on right pedal only	
	SRAM GX Shifter	
	SRAM GX Derailleur	
	Avid Brake Lever	
	Cable Driven Derailleur	
	11-36 Cassette	
TE65229	SRAM GS/GX Grip/Twist Shifter Package Shifter & brake on right pedal, left pedal mount available as a special	250,00
	SRAM GS Shifter	
	SRAM GX Derailleur	
	Avid Brake Lever	
	Cable Driven Derailleur	
	11-36 Cassette	
TE65067	SRAM Apex Road Shifting Package Shifter & brake mounted on right side of fork only	350,00
	SRAM Apex Shifter	
	SRAM GX Derailleur	
	Brake Lever Integrated into shifter	
	Cable Driven Derailleur	
	11-36 Cassette	
11 SPEED:		
TE65227	SRAM eTap Electronic Shifting Package Shifter on left pedal & brake on right pedal Not available with Triple chainring option (65104) Power Plate upgrade (65129) required"	1.750,00
	eTap Clic Shifter	
	SRAM Red eTap cassette Derailleur	
	SRAM Time Trial Brake Lever	
	SRAM Red eTap chainring Derailleur	
	11-36 Cassette	
BRAKING:		
TE65084	Rim Brakes, Tektro (see definitions)	
TE65082	Front disk brake Only Available with Wide Fork Option	120,00
TE65083	Rim Brakes, Torq Lite (see definitions)	450,00
CHAINRINGS:		
TE65104	Triple 48-36-26 Compact 110 Not available with eTap shifting (65227) SRM (65151)	
TE65059	Double 48-34 Compact 110	110,00
TE65106	Double 50-36 Compact 110	120,00
TRAINING SYSTEMS:		
TE65058	Omit SRM Training System	
TE65151	SRM Training System - Compact (Only Available with Double Chainrings)	3.900,00
FOOTREST OPTIONS: (see diagram below)		
TE60008	Standard Footrests, 3-3/4" Wide	
TE60009	Standard Footrests, 4-3/4" Wide	

	CODE	Omschrijving	Publieke prijs in €, 6% BTW inbegrepen
	TE60017	Connected Footrest 3-3/4" Wide with Bumper Not recommended for users under 5'8". Not available with the disc brake (65082) or wide fork (75026) options	50,00
	TE60018	Connected Footrest 4-3/4" Wide with Bumper Not recommended for users under 5'8". Not available with the disc brake (65082) or wide fork (75026) options	50,00
CRANKS: (see diagram below)			
CRANK STYLE:			
	TE65055	V Cranks	
CRANK LENGTH:			
	TE65036		160mm Long
	TE65037		165mm Long
	TE65038	170mm Long	
	TE65039	175mm Long	
	TE65146	Custom Crank Length	100,00
CRANK WIDTH:			
	TE65005	12" Wide	
	TE65007	13" Wide	
	TE65009	14" Wide	
	TE65011	15" Wide	
	TE65013	16" Wide	
	TE65015	17" Wide	
	TE65147	Custom Crank Width	100,00
UPHOLSTERY:			
MATERIAL:			
	TE30001	Nylon, Black Upholstery	
BACK UPHOLSTERY:			
	TE25064	Back Cushion	
SEAT UPHOLSTERY:			
	TE25063	Seat Cushion	
	TE25060	Hook and Loop Positioning Strap	
HAND PEDALS:			
SIZE:			
	TE65127	Medium (1-3/4") Ovalized Ergonomic	
	TE65113	Large (2-1/4") Ovalized Ergonomic	
ANGLE:			
	TE65171	15 Degrees	
	TE65172	25 Degrees Only Available with Apex (65067) or eTap (65227) shifting packages	
POWER PLATES:			
	TE65139	Omit Power Plates	
	TE65138	Power Plates Required with eTap shifting package (65227)	
FRAME COLORS:			
	55018		Wet Black 350,00
	55019	Black	Black Prism 350,00
	55017		Stealth Black 350,00
	55010		Midnight Blue 350,00
	55011	Blue	Royal Blue 350,00
	55012		Electric Blueberry 350,00
	55023	Grey	Stealth Grey 350,00
	55053		Agriculture Green 350,00
	55008	Green	Grasshopper Green 350,00
	55054		Neon Green 350,00

	CODE	Omschrijving	Publieke prijs in €, 6% BTW inbegrepen
	55006		Sunkissed Orange 350,00
	55005	Orange	Simple Orange 350,00
	55051		Neon Orange 350,00
	55015	Pink	Bubblegum Pink 350,00
	55057		Neon Pink 350,00
	55055	Purple	Lilac Purple 350,00
	55001	Red	Simple Red 350,00
	55002		Tangerine Red 350,00
	55003		Raspberry Red 350,00
	55050		Neon Red 350,00
	55020	Silver	Silver Metallic 350,00
	55016	White	White Snow 350,00
	55007	Yellow	Safety Yellow 350,00
	55052		Neon Yellow 350,00
	55066	Raw	Buffed Aluminum 350,00
FORK COLOR:			
Same as Frame Color (Not applicable with 55066, must select 55044)			
	55067	N/A	200,00
	55041	Black	Wet Black 200,00
	55042		Black Prism 200,00
	55040		Stealth Black 200,00
	55033	Blue	Midnight Blue 200,00
	55034		Royal Blue 200,00
	55035		Electric Blueberry 200,00
	55046	Grey	Stealth Grey 200,00
	55061	Green	Agriculture Green 200,00
	55031		Grasshopper Green 200,00
	55062		Neon Green 200,00
	55029	Orange	Sunkissed Orange 200,00
	55028		Simple Orange 200,00
	55059		Neon Orange 200,00
	55038	Pink	Bubblegum Pink 200,00
	55065		Neon Pink 200,00
	55063	Purple	Lilac Purple 200,00
	55024	Red	Simple Red 200,00
	55025		Tangerine Red 200,00
	55026		Raspberry Red 200,00
	55058		Neon Red 200,00
	55043	Silver	Silver Metallic 200,00
	55039	White	White Snow 200,00
	55030	Yellow	Safety Yellow 200,00
	55060		Neon Yellow 200,00
	55044	Raw	Buffed Aluminum 200,00
ACCESSORIES:			
STANDARD:			
	TE25061	Safety Flag	
	TE25067	Drafting Bumper	
	TE25070	Carbon Fiber Chainguard	
	TE25071	Adjustable Headrest	
	TE25038	Mirror; Chain Guard Mounted	
UPGRADABLE:			
SAFETY:			

	CODE	Omschrijving	Publieke prijs in €, 6% BTW inbegrepen
	TE25030	Extra Steering Dampener	
	TE25015	Light, Front (White)	
	TE25016	Light, Rear (Red)	
	TE25018	Helmet, Small/Medium	
	TE25019	Helmet, Medium/Large	
	TE25039	Leg Guard	
	TE25029	Tool and Tire Repair Kit	
	TE60015	Leg Support, Curved (Pair)	
POSITIONING:			
	TE25001	Click Strap, 8" Small	
	TE25002	Click Strap, 10" Medium	
	TE25003	Click Strap, 12" Large	
HYDRATION:			
	TE25021	Water Bottle and Cage	
	TE25020	Backpack Hydration System	
PERFORMANCE ENHANCEMENT:			
	TE25023	Crank Width Adaptor Set Not installed - (Gain +1.5" of crank width per set)	
	TE25022	Wheel Alignment Gauge	
	TE25024	Cordless Computer (tracks speed, distance, etc.)	
	TE25076	Garmin (tracks speed, distance, etc.)	
TRANSPORTATION:			
	TE25006	Crutch Holder	
	TE25007	Bike Rack, Hitch 1-1/4"	
	TE25008	Bike Rack, Hitch 2"	
BAGS:			
	TE25031	Wheel Bag	
	TE25032	Travel Bag	
TRAINING:			
	TE25033	Heartrate Monitor	
	TE25036	Indoor Stationary Roller	
	TE25037	Sportscrafter Mini Roller	
APPAREL:			
	TE25084	Lemon Yellow Top End T-Shirt	
	TE25085	Grey Top End T-Shirt	
	TE25086	Grey Top End Racer Back Tank Top	
Please select a size for your shirt selection:			
	TE25087	Extra Small	
	TE25088	Small	
	TE25089	Medium	
	TE25090	Large	
	TE25091	Medium	
	TE25092	XX Large	
	TE25093	XXX Large	
DECAL COLORS: (decals are reflective)			
FRAME DECAL COLOR:			
	TE55071	Black Frame Decals	
	TE55072	White Frame Decals	
FORK DECAL COLOR:			
	TE55073	Black Fork Decals	
	TE55074	White Fork Decals	

HANDCYCLE SPECIFICATIONS:	
Weight:	28 pounds
Turning Radius:	18 feet for 360 degree circle, 12 feet on a 3-point turn
Transfer Height:	5" from ground with 2" cushion
Frame Style:	Ovalized 7005 Aluminum tubing which is heat treated for increased strength and stiffness.
Fork:	Removable, Travel ready, with adjustable cranks
Propulsion /Shifting:	30 or 22 speeds with choice of components
Fork Style:	Adjustable length and height options
Back:	Carbon Fiber; adjustable angle from 50 to 70 degrees
Footrest:	Dual adjustable length, includes straps
Hand Pedals:	Ovalized ergonomic with frictionless bearings
Brake:	Caliper brakes standard, optional front disk brake
Wheelbase:	60" to 63" depending on wheel selection. See definitions
Camber:	3 or 0 degree available
Wheels:	20", 26" or 700c wheels with optional upgrades
Axles:	Threaded
Tires:	Standard tires provided. Exact brand varies based on availability
Positioning Straps:	Hook and loop positioning strap included
Cushions:	Seat and back cushion included
Safety / Visibility:	Flag and reflectors included
Chain Guard:	Carbon Fiber chain guard included

EXCLUSIONS AND DEFINITIONS:	
Top End Force RX Frame and Fork:	The Top End Force RX handcycle is made of 7005 series aluminum. An adjustable carbon fiber backrest and an adjustable bolt-on headrest are included. With a standard configuration the Force RX is 30 speeds and with eTap the bike is 22 speeds.
Assembled vs Unassembled and shipping:	If TE10046 (Unassembled) is selected the bike will require about 60 minutes of assembly upon delivery (instructions provided). If TE10047 (Assembled) is selected the bike will be fully assembled upon delivery. Shipping rate is \$150 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
Wheel Sets:	Wheel sets include a front wheel, two rear wheels with axles, tires and tubes. Wheel sets are not able to be mixed and matched or broken out into individual wheel selections.
Tires:	Must choose the same size tires as wheels. Clincher wheels must have clincher tires and tubular wheels must have tubular tires.
Fork Style:	Fork adjustment range is about 19" to 21.5" with a 26" front wheel and 20" to 22.5" with a 700c front wheels
Connected Footrests:	The connected footrests are not recommended for athletes under 5'8" (1.77m) and are not available with the disc brake (65082) or wide fork (75026) options.
eTap Electronic Shifting:	Triple chainring combination not available, must select double chainring options 65059 or 65106. Power Plate upgrade required.
Power Plates:	Required with eTap shifting package (65227)
Omit Wheels:	When an omit wheels option is selected the following will NOT be included. Front or rear wheels, tires, tubes, axles. When a shifting package is selected the cassette from that package will be included.
Wheelbase:	Wheelbase is measured from the middle of the rear axle to the middle of the front hub. When Standard length is chosen wheelbase will be 60" with a 26" front wheel and 61" with a 700c front wheel. When extended length is chosen wheelbase will be 61" with 26" front wheel and 63" with a 700c front wheel.
Braking Options:	The Tektro and Torq Lite braking packages include a caliper and cables for the hand brake lever included in the shifting package and flip lever parking brake assembly. The disc brake braking package includes a disc brake rotor, caliper, cables and a Tektro flip lever parking brake assembly.
Omit Shifting:	When omit shifting option is selected the shifter, chainring derailleur, cassette derailleur, cassette, friction shifter for brake caliper, chain and shifter cables will NOT be included. The handpedals and fork will be made at accept trigger shifter components mounted on the right handpedal. A hand brake lever, caliper and cables as well as a flip-lever parking brake will be included. The brand/style is determined by the braking option selected.
25 Degree Handpedal Angle:"	Only available with the Apex fork mounted (65067) or the eTap (65227) shifting packages due to interference.
SRM Training System:	Only available with double chainring options. Not available with Triple chainring option 65104.
Axles:	
Tires:	
Positioning Straps:	
Cushions:	
Safety / Visibility:	
Chain Guard:	

Warranty Info:

Invacare warrants the frames and forks when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Invacare or a dealer; with a copy of the seller's invoice required for coverage under this warranty. Invacare warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Invacare or a dealer; with a copy of the seller's invoice required for coverage under this warranty

Note: All Special requests are subject to approval, additional charges, and additional lead time

Special Request 1)

Special Request 2)

Special Request 3)

Special Request 4)

Special Request 5)