


Customer :  
 Customer No. :  
 Phone :

Contact :  
 Commission :  
**Serial No. :**



Item	Qty	Part-No.	Part-Designation
-	1	SP1428069	BATTERY COVER BLUE CPL.
-	1	SP1421846	BATTERY COVER SILVER CPL.
-	1	SP1421845	BATTERY COVER GREEN CPL.
1	1		INVACARE LOGO
2	1		LOOPED VELCRO STRIP 25X250 MM, BLACK
-	1		HOOKED VELCRO STRIP 25X250 MM, BLACK
3	1		BATTERY COVERING "SILVER"
3	1		BATTERY COVERING "GREEN"
3	1		BATTERY COVER "BLUE"
-	1		BATTERY CABLE DIAGRAM STICKER
4	2	P50016	BATTERY 12V / 70 AH
-	1	1421655	BATTERY STRAP
-	4	P52996	BATTERY TERMINAL COVER
-	3	F384057	BATTERY BOX LINING, BLACK
-	1	<b>SP1423696</b>	<b>CHARGER 24V/8A W. CONTINENTAL MAINS CABLE, CPL.</b>
5	1	P51981	BATTERY CHARGER 8 A WITH XLR CONNECTION
6	1	P21079	MAINS CABLE GCB-0401 CONTINENTAL; BLACK
-	1	<b>SP1423697</b>	<b>CHARGER 24V/8A W. UK MAINS CABLE, CPL.</b>
5	1	P51981	BATTERY CHARGER 8 A WITH XLR CONNECTION
-	1	P21080	MAINS CABLE GCB-0403 UK; BLACK

### Storm 3

