

Spare part list OCEAN VIP

No.	Qty	Item	Art.-Nr.
1	1	Brake cable assembled (all cables for tilting-mechanism + divider)	1470685
2	1	Clip footrest	16333
3	1	Clip frame	1470686
4	2	Heel loop (pair)	16334
5	1	5" castor with directional lock	1470687
6	1	5" castor	16259
7	1	Threaded plate	1470688
8	1	Footrest plate	1470689
9	1	Seat	1470690
10	1	Headrest bracket, set, new design	1470971
11	1	Fastening set, backrest (four screws and washers)	16434
12	1	Backrest cover	1470692
13	1	Armrest assembly, including holder, pre-assembled, right, including screws	1470693
14	1	Armrest assembly, including holder, pre-assembled, left including screws	1470694
15	1	Footrest complete, left (tubes, footrest plate, insert, heel loop)	1470695
16	1	Footrest complete, right (tubes, footrest plate, insert, heel loop)	1470696
17	1	Armrest pad, including screws	16472
18	1	Plugs for sideparts (rear), set (8 pcs)	1470698
19	1	Adapter for backrest, triangles including screws	1470699
20	1	Headrest complete, assembled w/o bracket	1470700
21	1	Gasspring	1470701
22	1	Plugs for sideparts (front), set (6 pcs)	1470702
23	1	Screws for basic tube, set	1470704
24	1	Release mechanism for gasspring	1470703
25	1	Brake lever	1470732
26	1	Headrest complete, including bracket	1470790
27	1	Footrestholder, new design	1470929
28	1	Insert for headrest bracket	1514509