

.....

Spare Part List - INVACARE® Storm³



Table of Contents - INVACARE Storm³

.....	1
Warnings/ Information	5
Storm ³	7
Dataplate	7
Chassis	8
Crash Kit	9
Shrouds	11
Wheel Suspension/ Wheels	12
Front Wheels	14
Driving Wheels	15
Rear Wheels True Track	17
Driving unit	18
Brake	20
Decoupling	21
Electric Support	22
Lighting	23
Front Lights	23
Rear Lights	24
Electronic	25
Cables	26
ACS-Buscables	26
CLAM-Cables (LIGHTING/ACTUATOR MODULE COMBINED)	26
Battery Cables	26
Legrests, sideways mounted	27
Footrest receptacle	27
Footrest receptacle complete, manual 8cm	28
Footrest receptacle cpl., electric 8cm	29
Footrest receptacle, complete, manual 5cm	30
Footrest receptacle, cpl., electric 5cm	31
Hanger	32
Vari F+A Detailed	33
Release unit assembly kit	34
Upper Vari F, right	35
Upper Vari F, left	36
Release unit assembly kit f. legrest w. man. height adjustment	37
Upper Vari A complete, right	38
Upper Vari A complete, left	39
ADM + ADP Detailed	40
Release unit assembly kit f. legrest 80° w. man. height adjustment	41
Complete upper ADM, right	42
Assembly kit, ADM hangers	43
Assembly kit, ADM sliding bush	44
Assembly kit, Trip lever + lock lever	45
Assembly kit, ADM adjustment	46
Complete upper ADM, left	47
Assembly kit, ADM hangers	48
Assembly kit, ADM sliding bush	49
Assembly kit, Trip lever + lock lever	50
Assembly kit, ADM adjustment	51
Upper ADP, right	52
Assembly kit, ADM hangers	53
Assembly kit, ADM sliding bush	54

Assembly kit, Trip lever + lock lever.....	55
Assembly kit ADP, actuator bracket.....	56
Assembly kit ADP, hoist motor.....	57
Upper ADP, left.....	58
Assembly kit, ADM hangers.....	59
Assembly kit, ADM sliding bush.....	60
Assembly kit, Trip lever + lock lever.....	61
Assembly kit ADP, actuator bracket.....	62
Assembly kit ADP, hoist motor.....	63
Lower Footrests.....	64
Lower Legrests Detailed.....	65
Calf plate holder, right.....	66
Calf plate holder, left.....	67
Back Unit.....	68
Backrest frame Standard.....	68
Backrest frame CONTOUR.....	69
Backrest frame FLEX.....	70
Backrest adjustment STANDARD.....	71
Backrest adjustment FLEX/CONTUR man.....	72
Backrest adjustment FLEX/CONTUR el.....	73
Backrest adjustment ASSEMBLY KIT, cpl.....	75
Retaining belts.....	76
Head rests.....	77
Backrest Upholstery.....	79
Side Panels.....	81
Euro Side Panels.....	81
Standard Side Panels.....	82
Side Panels for Standard seatheight.....	83
Side part plates complete with side part adjustment.....	85
Side part receptacles complete.....	86
Side part receptacles single.....	87
Side Panels for Arne Jensen-Seats.....	88
Sliding Armrests.....	90
Lifter.....	91
Seat Height Adjustment Manual.....	92
Seat Angle Adjustment Manual.....	93
Seat Height Adjustment electrical.....	94
Seat Angle Adjustment electrical.....	95
Control Unit/ Joystickbox.....	96
Attendant Control.....	98
Central Steering.....	99
Chin Control Dash Board.....	100
Chin Control Assembly Kit.....	101
Control Unit support.....	102
Standard/ swing away.....	102
For Attendant Conroll.....	103
Control Unit support tabular.....	104
Desk/ Tray control.....	105
Tray Control.....	106
Options.....	107
Kerb Climber.....	108
Documentation.....	109
User Manuals Storm ³	109
User Manuals REM 24 SC -> 03/2004.....	110
User Manuals REM 24 SD.....	111

User Manuals SHARK.....	112
User Manuals miscellaneous	113

Warnings/ Information

⚠ **Caution: Physical injury and damage to property! Loss of CE conformity and warranty due to modifications and mounting of privately purchased articles without authority!**

Alterations to the vehicle that are not approved by Invacare® and/or do not belong to the standard or optional equipment can

- downgrade vehicle handling to such an extent that the vehicle is no longer controllable,
- restrict functionality to such an extent that this may cause failure of electric/electronic components or become a fire hazard.

Modifications that have not been approved and the mounting of privately purchased articles will lead to the loss of CE conformity, Invacare® will not be liable in the event of physical injury or damage to property.

Expiration of vehicle warranty.

Maintain full efficiency and the safe use of your vehicle, avoid physical injury, damage to property or the loss of CE conformity and warranty by:

- only using original Invacare® spare parts,
- only using installation material of the same quality as the original material,
- never opening electric and electronic components (not even after warranty has expired!),
- never modifying electric, electronic or mechanical components
- never mounting parts obtained from external manufacturers without written approval from Invacare®.

✓ **Information: Definitions.**

What parts are looked upon by Invacare® as being from external manufacturers and as electric/electronic components?

Among other things, parts from external manufacturers that are subject to approval are:

- All electric and electronic components that are connected to the power supply system of the vehicle, such as seat heating, remotes, signal horns, radios.
- All mechanical mounted parts that change vehicle handling, such as trailed units and seat systems.
- In general mounted parts are not permitted that do not conform with the valid road traffic regulations and may thus also not be mounted on Invacare® vehicles.

Electric and electronic components must not be opened, e.g.:

- Electronics
- Drive motors
- Remotes
- Actuators for electric adjustments, e.g. adjustment of the back, legrests and seat angle.

⚠ **Caution: Risk of injury when assembling and dismantling the vehicle!**

Before commencing any work, ensure that the vehicle is secure. Please read and follow the valid industrial safety regulations. Carry out maintenance and service work according to the service instructions (available separately) in order to avoid physical injury and damage to property. A tool list, illustrated step-by-step instructions and details of tightening torques, will assist you and show how the vehicle remains, or will remain safe and fully serviceable.

✓ **Information: Invacare® addresses and service numbers.**

Addresses and service numbers of Invacare® branches can be found in the operating instructions, in the service instructions and on the Internet.

✓ **Information: Meaning of symbols.**

⚠ **Caution** ✓ **Information**



LH - left hand



RH - right hand



black



grey



Vinyl grey



Dartex grey



blue



blue birdie



blue dots



green birdie



green dots


Dataplate



- ✓ You are often only able to find the suitable spare part by using the details on the type plate. The type plate is to be found directly on the chassis.

Please quote the serial number, model and year for every order. This helps to prevent misunderstandings.

Typenschild.dwg



item	description
1	MODEL DESIGNATION
2	SERIAL NUMBER
3	PRODUCTION DATE