



## Remote Holder

G50\_Fahrpulhalter.dwg



item	ordercode	description	quantity
1	SPS183040	REMOTE HOLDER STANDARD	1
2	L31001	HAND KNOB M6 X 10	1
3	SP1421796	SWING-OUT ACS REMOTE HOLDER LH	1
3	SP1421795	SWING-OUT ACS REMOTE HOLDER RH	1
4	SP1422937	ASSEMBLY KIT F. STD. CARER SEAT CONTROL LH .	1
4	SP1422936	ASSEMBLY KIT RH F. STD. CARER SEAT CONTROL	1
5	SP1422819	ASSEMBLY KIT F. SWITCHABLE CARER RECARO CONTROL	1
6	SP1422818	ASSEMBLY KIT F. SWITCHABLE CARER CONTROL CPL.	1
7	1440631	TUBE CLAMP DIA 22 MM; BLACK	1
8	Z183264.99	BUSH RD. 22 X 1,5 X 40	4
9	Z12008A.P40	GRIPPING CLIP W. HEXAGON HOLE	2
10	L31001	HAND KNOB M6 X 10	1
11	F383065-9005	HOLDER FOR CARER RECARO CONTROL; RAL 9005, DEEP BLACK	1
12	SP1440289	ATTENDANT JOYSTICK BOX SUPPORT REM 24SD WITH SCREWS, LEFT	1
12	SP1440288	ATTENDANT JOYSTICK BOX SUPPORT REM 24SD WITH SCREWS, RIGHT	1
-	SP1422814	ASSEMBLY KIT F. SWITCHABLE CARER CONTROL CPL., RECARO LH	1
-	SP1422813	ASSEMBLY KIT F. SWITCHABLE CARER CONTROL CPL., RECARO RH	1
13	SP1434921	JOYSTICK BOX SUPPORT SWING-OUT REM 24SD AJ, WITH SCREWS	1
14	SP1436542	REMOTE HOLDER REM 24SD	1

## Remote Holder for a accompany person



item	ordercode	description	quantity
1	SP1422937	ASSEMBLY KIT F. STD. CARER SEAT CONTROL LH .	1
1	SP1422936	ASSEMBLY KIT RH F. STD. CARER SEAT CONTROL	1
2	SP1422818	ASSEMBLY KIT F. SWITCHABLE CARER CONTROL CPL.	1
3	SP1422814	ASSEMBLY KIT F. SWITCHABLE CARER CONTROL CPL., RECARO LH	1
3	SP1422813	ASSEMBLY KIT F. SWITCHABLE CARER CONTROL CPL., RECARO RH	1
4	SP1422819	ASSEMBLY KIT F. SWITCHABLE CARER RECARO CONTROL	1
5	L31001	HAND KNOB M6 X 10	1
6	Z12008A.P40	GRIPPING CLIP W. HEXAGON HOLE	2
7	F383065-9005	HOLDER FOR CARER RECARO CONTROL; RAL 9005, DEEP BLACK	1
8	SP1440289	ATTENDANT JOYSTICK BOX SUPPORT REM 24SD WITH SCREWS, LEFT	1
8	SP1440288	ATTENDANT JOYSTICK BOX SUPPORT REM 24SD WITH SCREWS, RIGHT	1